

# Feature extraction

Robert Haase

With material from


Johannes Müller, PoL TU Dresden

Marcelo Zoccoler, PoL, TU Dresden


Benoit Lombardot, Scientific Computing Facility, MPI CBG

December 2022

- Feature extraction is a *late* processing step in image analysis.
- It can be used for images, or segmented/labelled images



Feature  
Extraction



Acquisition






Denoising

Background subtraction

Segmentation




Labeling

Feature  
Extraction




- A *feature* is a countable or measurable property of an image or object.
- Goal of feature extraction is finding a minimal set of features to describe an object well enough to differentiate it from other objects.
- **Intensity based**
  - Mean intensity
  - Standard deviation
  - Total intensity
  - Textures
- **Shape based /spatial**
  - Area / Volume
  - Roundness
  - Solidity
  - Circularity / Sphericity
  - Elongation
  - Centroid
  - Bounding box
- **Spatio-temporal**
  - Displacement,
  - Speed,
  - Acceleration
- **Topological**
  - Number of neighbors
- **Others**
  - Overlap
  - Colocalization
- **Mixed features**
  - Center of mass
  - Local minima / maxima
  - Distance to neighbors
  - Average intensity in neighborhood


- Min / max
- Median
- Mean
- Mode
- Variance
- Standard deviation




- Can be derived from pixel values
- Don't take spatial relationship of pixels into account
- See also:
  - descriptive statistics
  - histogram



Count: 65024  
Mean: 103.301  
StdDev: 57.991  
Min: 29  
Max: 248  
Mode: 53 (1663)



Count: 783  
Mean: 141.308  
StdDev: 61.876  
Min: 44  
Max: 243  
Mode: 236 (9)



Count: 1056  
Mean: 49.016  
StdDev: 12.685  
Min: 34  
Max: 122  
Mode: 45 (120)

- Relative position in an image weighted by pixel intensities


- $x, y$  ... pixel coordinates
- $w$  ... image width
- $h$  ... image height
- $\mu$  ... mean intensity
- $g_{x,y}$  ... pixel grey value
- $x_m, y_m$  ... center of mass coordinates

$$\mu = \frac{1}{wh} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} g_{x,y}$$

$$x_m = \frac{1}{wh\mu} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} x g_{x,y}$$

$$y_m = \frac{1}{wh\mu} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} y g_{x,y}$$

“sum intensity”  
“total intensity”



$$x_m = 1/7 (1 \cdot 0 + 1 \cdot 1 + 2 \cdot 2 + 2 \cdot 1 + 1 \cdot 2) = 1.3$$

$$y_m = 1/7 (1 \cdot 2 + 1 \cdot 2 + 2 \cdot 3 + 2 \cdot 2 + 1 \cdot 3) = 2.4$$

- Relative position in an image weighted by pixel intensities
- Special case of center of mass for binary images
  - $x, y$  ... pixel coordinates
  - $w$  ... image width
  - $h$  ... image height
  - $\mu$  ... mean intensity
  - $g_{x,y}$  ... pixel grey value, integer in range [0;1]
  - $x_m, y_m$  ... center of mass coordinates

$$\mu = \frac{1}{wh} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} g_{x,y}$$

$$x_m = \frac{1}{wh\mu} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} x g_{x,y}$$

$$y_m = \frac{1}{wh\mu} \sum_{y=0}^{h-1} \sum_{x=0}^{w-1} y g_{x,y}$$



Number of white pixels

	0	1	2	3	4	x
0	0	0	0	0	0	
1	0	0	0	0	0	
2	1	1	1	0	0	
3	0	1	1	0	0	
4	0	0	0	0	0	
y						




$$x_m = 1/5 (1 \cdot 0 + 1 \cdot 1 + 1 \cdot 2 + 1 \cdot 1 + 1 \cdot 2) = 1.2$$

$$y_m = 1/5 (1 \cdot 2 + 1 \cdot 2 + 1 \cdot 3 + 1 \cdot 2 + 1 \cdot 3) = 2.4$$

- Length of the outline around an object
- Depends on the actual implementation




- For every object, find the optimal ellipse simplifying the object.
- Major axis ... long diameter
- Minor axis ... short diameter
- Major and minor axis are perpendicular to each other






- The aspect ratio describes the elongation of an object.


$$AR = \text{major} / \text{minor}$$




- By removing all concave corners of an object, we retrieve its **convex hull**.




- By removing all concave corners of an object, we retrieve its **convex hull**.



- By removing all concave corners of an object, we retrieve its **convex hull**.



- By removing all concave corners of an object, we retrieve its **convex hull**.



$$solidity = \frac{A}{A_{convexHull}}$$

# Roundness and circularity

- The definition of a circle leads us to measurements of circularity and roundness.
- In case you use these measures, define them correctly. They are not standardized!

Diameter


$d$

Circumference

$$C = \pi d$$

Area

$$A = \frac{\pi d^2}{4}$$



$$roundness = \frac{4 * A}{\pi major^2}$$


$$circularity = \frac{4\pi * A}{perimeter^2}$$

Roundness = 1  
Circularity = 1

Roundness  $\approx$  1  
Circularity  $\approx$  1

Roundness < 1  
Circularity < 1

# Roundness and circularity



$$\text{roundness} = \frac{4 * A}{\pi \text{major}^2}$$

$$\text{circularity} = \frac{4\pi * A}{\text{perimeter}^2}$$

$$\text{solidity} = \frac{A}{A_{\text{convexHull}}}$$

- In Python: `from skimage import measure`

<https://scikit-image.org/docs/stable/api/skimage.measure.html>

`skimage.measure.blur_effect` (image[, h\_size, ...]) Compute a metric that indicates the strength of blur in an image (0 for no blur, 1 for maximal blur).

`skimage.measure.euler_number` (image[, ...]) Calculate the Euler characteristic in binary image.

`skimage.measure.find_contours` (image[, ...]) Find iso-valued contours in a 2D array for a given level value.

`skimage.measure.grid_points_in_poly` (shape, verts) Test whether points on a specified grid are inside a polygon.

`skimage.measure.inertia_tensor` (image[, mu]) Compute the inertia tensor of the input image.

`skimage.measure.inertia_tensor_eigvals` (image) Compute the eigenvalues of the inertia tensor of the image.

`skimage.measure.label` (label\_image[, ...]) Label connected regions of an integer array.

`skimage.measure.regionprops` (label\_image[, ...]) Measure properties of labeled image regions.

`skimage.measure.regionprops_table` (label\_image) Compute image properties and return them as a pandas-compatible table.

**area** : int

Number of pixels of the region.

**area\_bbox** : int

Number of pixels of bounding box.

**area\_convex** : int

Number of pixels of convex hull image, which is the smallest convex polygon that e

**area\_filled** : int

Number of pixels of the region will all the holes filled in. Describes the area of the i

**axis\_major\_length** : float

The length of the major axis of the ellipse that has the same normalized second ce  
the region.

**axis\_minor\_length** : float

The length of the minor axis of the ellipse that has the same normalized second ce  
the region.









# Exploring neighborhood relationships between cells

- Study how many neighbors objects have.
- How likely is it that an object with 3 neighbors is a cell?


```
num_neighbors_map = cle.touching_neighbor_count_map(objects)  
  
cle.imshow(num_neighbors_map,  
           color_map='jet',  
           colorbar=True)
```

```
cle.imshow(num_neighbors_map,  
           min_display_intensity=0,  
           max_display_intensity=6,  
           color_map='jet',  
           colorbar=True)
```




- Filter out objects which have an unreasonable number of neighbors


```
cells = cle.exclude_labels_with_map_values_out_of_range(num_neighbors_map, objects, minimum_value_range=4)  
cle.imshow(cells, labels=True)
```



neighbor count




objects



cells

- Filter labeled objects using Measure Labels and Label Filters in Napari.



# SimpleITK

- Recommended for 3D-measurements, based on the SimpleITK-project

The screenshot displays the napari software interface with three main windows and a properties panel.

**Left Window (napari):** Shows a 3D visualization of segmented cells. The layer controls on the left include settings for label, opacity (0.7), brush size (10), blending (translucent), and color mode (auto). The layer list on the bottom left shows 'connected\_compo...', 'threshold\_otsu res...', and 'median\_filter result'.

**Middle Window (napari):** Shows a 2D view of the same data. The layer controls on the left show 'Threshold (Otsu et al 1979, n-SimpleITK)' with 'image' set to 'median\_filter result'. The layer list on the bottom left shows 'connected\_compo...', 'threshold\_otsu res...', and 'median\_filter result'.

**Right Window (napari):** Shows a 2D view of the segmented cells. The layer controls on the left show 'Threshold (Otsu et al 1979, n-SimpleITK)' with 'image' set to 'blobs (data)'. The layer list on the bottom left shows 'connected\_compo...', 'threshold\_otsu res...', and 'median\_filter result'.

**Properties Panel (Properties of Result of connected\_component\_labeling):** A table showing the results of the connected component labeling process.

label	maximum	mean	median	minimum
1	232.0	190.8545034642...	200.0	128.0
2	224.0	179.2864864864...	184.0	128.0
3	248.0	205.6170212765...	208.0	128.0
4	248.0	217.3271889400...	232.0	128.0
5	248.0	212.1425576519...	224.0	128.0
6	248.0	204.2947368421...	216.0	128.0
7	200.0	161.4814814814...	168.0	128.0

- Many Napari plugins for feature extraction can also be called from Python

```
statistics = label_statistics(blobs, labels,  
                             intensity=True,  
                             size=True,  
                             shape=True,  
                             perimeter=True,  
                             position=True,  
                             moments=True)
```


```
df = pd.DataFrame(statistics)  
df
```

	label	maximum	mean	median	minimum	sigma	sum	variance	bbox_0	bbox_1
0	1	224.0	137.526132	136.0	112.0	13.360739	157880.0	178.509343	0	0
1	2	232.0	193.014354	200.0	128.0	28.559077	80680.0	815.620897	11	0
2	3	224.0	179.846995	184.0	128.0	21.328889	32912.0	454.921516	53	0
3	4	248.0	207.082171	216.0	120.0	27.772832	133568.0	771.330194	95	0
4	5	248.0	223.146402	232.0	128.0	30.246515	89928.0	914.851647	144	0
5	6	248.0	214.906725	224.0	128.0	26.386796	99072.0	696.263020	238	0
6	7	248.0	211.565891	224.0	136.0	30.197236	54584.0	911.873073	189	7
7	8	200.0	166.171429	168.0	136.0	16.466894	11632.0	271.158592	133	17




# Exploring features in Napari

- Select table rows and view corresponding object in 2D/3D space



# Exploring features in Napari

- Double-click on table column to retrieve a parametric map image



# Exercise: Quantitative measurements

- Use the given feature extraction notebook to apply some basic statistics to measurements

```
# analyse objects
properties = measure.regionprops(label_image, intensity_image=image)
```

```
statistics = {
    'area':      [p.area          for p in properties],
    'mean':      [p.mean_intensity for p in properties],
    'major_axis': [p.major_axis_length for p in properties]
}
```

```
dataframe = pd.DataFrame(statistics)
dataframe
```

	area	mean	major_axis	aspect_ratio
0	429	191.440559	34.779230	2.088249
1	183	179.846995	20.950530	1.782168
2	658	205.604863	30.198484	1.067734
3	433	217.515012	24.508791	1.061942
4	472	213.033898	31.084766	1.579415
...	...	...	...	...
57	213	184.525822	18.753879	1.296143
58	79	184.810127	18.287489	3.173540
59	88	182.727273	21.673692	4.021193
60	52	189.538462	14.335104	2.839825
61	48	173.833333	16.925660	4.417297

62 rows × 4 columns

- How many objects are in it?

- How large is the largest object?

- What are mean and standard deviation of the intensity in the image?

- What are mean and standard deviation of the area of the segmented objects?



# Exercise: Parametric maps

- Produce a parametric map representing 'elongation' in Napari.


The image displays three sequential screenshots of the Napari software interface, illustrating the steps to create a parametric map for elongation.

**Left Screenshot:** Shows the initial state with a grayscale image of cells. The 'layer controls' panel on the left is visible, showing the 'label' set to 1 and 'opacity' at 0.7. The 'Measurement' panel on the right shows the 'blobs' layer selected, with 'shape' checked under the 'Operation' dropdown.

**Middle Screenshot:** Shows the 'Properties of Result of Connected component labeling' dialog box. The 'Copy to clipboard' button is highlighted. The 'label' column is visible, showing values 1, 2, 3, and 4.

**Right Screenshot:** Shows the final state with a color-coded parametric map. The 'layer controls' panel on the left shows the 'opacity' set to 1.0. The 'Measurement' panel on the right shows the 'blobs' layer selected, with 'shape' checked under the 'Operation' dropdown. The 'Properties of Result of Connected component labeling' dialog box is also visible, showing the 'Copy to clipboard' button and the 'label' column.

use <1> for activate the label eraser, use <2> for activate the paint brush, use <3> for activate the fill bucket, use



```
filtered = filters.median(image)  
  
filtered = filters.gaussian(image, sigma=5)
```

Filtering the image reduces pixel noise

```
bg_subtracted = morphology.white_tophat(image, footprint=footprint)
```

Top-hat filtering removes the background

Thresholding binarizes the image

```
threshold = filters.threshold_otsu(image)
```

Connected-components analysis groups pixels to objects

```
labels = measure.label(binary)
```

Feature extraction allows descriptive statistics

```
measurements = measure.regionprops_table(labels, properties=properties)
```